EZ-Fit® Translucent

Multi-Tiered, Translucent
Passive Fiber Post

- Size 0 Posts
- Size 1 Posts
- Size 2 Posts

Essential Dental Systems, Inc.
89 Leuning Street, S. Hackensack, NJ 07606, U.S.A.

Copyright © 2014 Essential Dental Systems, Inc. All rights reserved.

Sterilization Procedure:

Posts, Drills and Wrenches are non-sterile.

Prior to use, bag materials and sterilize with a fractionated vacuum autoclave or gravity displacement autoclave for 20 minutes at 121°C. When reprocessing drills, remove debris and dry before sterilization. Remove debris with an ultrasonic cleaning unit for 5 minutes or a thermal disinfectant unit (washer-disinfectant) at 80°C for 10 minutes. For complete instructions see website.

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ENGLISH: Post selection is best done by placing a post in a hemostat and holding it over an undistorted x-ray. If there is a minimum of approximately 0.5mm of lateral root structure at the most apical placement of the post, there is sufficient root structure for placement. Make sure to use a large enough post size (i.e. diameter) to withstand the occlusal forces. Note that the EZ-Fit Post is designed with a second tier for reduction of apical stress. The second tier provides increased retention at the critical coronal level (where the canal flares), stabilizing the post for a long lasting restoration.

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1. ENGLISH: Use a Peeso or Essential Gates Glidden drill to remove gutta percha, and for preliminary sizing of the canal (EDS' Gates Glidden Drills are recommended). Then use the primary reamer.

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2. ENGLISH: Use the primary reamer to prepare the full length of the post hole (Cut Wet).

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3. ENGLISH: Use the secondary drill to prepare the seat for the second tier. (The second tier must be fully seated. If not, you reduce the post's retention and increase the chance of fracture under function.)

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


4. ENGLISH: Trial insert the post to determine fit. Note second tier must fully seat.

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
5. ENGLISH: Post adjustment: cut apical end, if necessary, to ensure seating of the second tier. Note: apical adjustment can be done after trial insertion.

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6. ENGLISH: Note: For optimal retention EDS recommends the use of EDTA-Plus (EDS SKU 770-02) followed by the use of EDS' Bonding Agent (EDS SKU 900-00) and Envy Cement EDS SKU 960-00).

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
7. ENGLISH: Dry the canal. Place cement in the canal. **Note fully seated second tier (A).** (Envy[®] Cement from EDS is recommended). Follow with core formation and preparation (Ti-Core[®] composite. Core build-up material from EDS is recommended).

As the post is highly translucent, you may light cure through the post, down the canal. If doing so, wipe away any excess cement before activating the light and follow manufacturer's instruction. Place curing light probe directly over the post, applying gentle pressure with the curing light tip. Light-cure for a minimum of 60 seconds (depending on volume and depth).

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Notes: Keep posts out of sunlight. In all EZ-Fit post restorations, a minimum of a 1.5mm ferrule around the tooth is needed for the final restoration.

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FOR ADDITIONAL PRODUCT INFORMATION AND ORDERING PLEASE VISIT: WWW.EDSDENTAL.COM

WARNING: This product can expose you to the chemical Bisphenol A, which is known to the State of California to cause reproductive harm. For more information go to www.P65Warnings.ca.gov.

EDS ESSENTIAL[®] DENTAL SYSTEMS

89 Leuning Street, S. Hackensack, NJ 07606, U.S.A.

1-800-22-FLEXI

Tel: (201) 487-9090

Fax: (201) 487-5120

www.edsdental.com

EZFITFRANS-1

ENDODONTIC FIBER POST CLEANING

While chemical sterilants such as peracetic acid, 2% glutaraldehyde based formulations, and/or 6% hydrogen peroxide may be used for cleaning endodontic fiber posts, *the majority of fiber endodontic post manufacturers recommend the use of disinfectants such as 70% ethyl or isopropyl alcohol*^{1,2,3,4,8}. This has become a dental industry standard for the cleaning of all fiber post systems before placement. The following protocol is recommended by most dental industry's top manufacturers including 3M, Brassler, as well as Dentsply Sirona^{5,6,7}:

- Disinfect the post with alcohol and dry it with air free of water and oil.
- Wipe the post with alcohol.
- Allow the post to air dry for 30 seconds before use.

While not necessary, Endodontic glass fiber posts may also be sterilized using conventional methods according to a study by *Yagci et al^a without any deleterious structural effects to the post*. In this study, conventional steam autoclave sterilization was performed on a large number of glass fiber posts and then the bond strength, flexural strength as well as the elastic modulus of each was evaluated. Additionally, the microstructures of the tested posts were evaluated by SEM and energy-dispersive X-ray spectroscopy (EDS) for damage resulting from the autoclave steam sterilization process. *One-way ANOVA indicated no statistically significant differences* among the groups means in terms of bond strength ($p > 0.05$), flexural strength ($p > 0.05$), or elastic modulus ($p > 0.05$). *The results of this study show that glass-fiber posts can be safely and effectively sterilized either by autoclave when necessary, without any negative effect on bond strength, flexural strength, or elastic modulus*⁴.

While studies have shown that glass fiber posts are not negatively affected during the steam autoclave sterilization process^{2,3,4,8}, cleaning of glass fiber posts with alcohol ($\geq 70\%$) is adequate for proper disinfection and recommended as per industry standard protocols^{5,6,7}. The cleaning, disinfection, and sterilization of glass fiber posts is ultimately *left to the discretion of the clinician* as dental industry standard protocol and/or steam autoclave sterilization procedures are both effective/safe. The disinfection of glass fiber posts using alcohol has long been an industry standard with an extensive history of safe clinical use and efficacy.

Citations

1. Bateman G, Ricketts DNJ, Saunders WP. Fiber-based post systems: a review. Br Dent J 2003;195:43–48.
2. Farrugia C, Cassar G, Valdramidis V, Camilleri J. Effect of sterilization techniques prior to antimicrobial testing on physical properties of dental restorative materials. J Dent 2015;43:703–714.
3. Kotha SB, Ramakrishnaiah R, Devang DD, Celur SL, Qasim S, Matinlinna JP. Effect of disinfection and sterilization on the tensile strength, surface roughness, and wettability of elastomers. J Investig Clin Dent 2016;8: doi:10.1111/jicd.12244.
4. Yagci F, Ustun Y, Zortuk M, Agirnasligil M. Effect of Sterilization on Bond Strength and Mechanical Properties of Fiber Posts. J Adhes Dent. 2019;21(2):143-148. doi: 10.3290/j.jad.a42325. PMID: 30949627.
5. <https://multimedia.3m.com/mws/media/3858070/3m-relyx-fiber-post-cementing-endodontic-post.pdf>.
6. <http://brasslerusa.com/wp-content/files/IFU-0026%20ES%20Fiber%20Post%20System%20IFU%20REV%20A.pdf>
7. <https://www.dentsplysirona.com/content/dam/master/product-procedure-brand-categories/endodontics/product-categories/ifu/END-IFU-Endodontic-Sterilization-Procedure-EN-1402.pdf>
8. Patel SS, Sethuraman R. A randomized controlled twelve month clinical study on the evaluation of success rate of endodontically treated teeth restored with metal poly-fiber posts and dentin posts. J Indian Prosthodont Soc. 2022 Jan-Mar;22(1):38-47. doi: 10.4103/jips.jips_134_21. PMID: 36510946; PMCID: PMC8884352.