ENGLISH: SafeSiders® HF Instructions for use.

NOTE: The SafeSiders® HF instruments have been designed for use in Rotation and for optimal use in the Endo-CSV™ Endodontic Handpiece. They should never be used in Asymmetric Oscillation/Reciprocation motors or handpieces.

For Use with The Endo-CSV Endodontic Handpiece:

Indications for use: For the removal of dentin and shaping of the root canal.

Contraindications: None Known.

Warnings:

- Instruments are packaged non-sterile, and must be sterilized before use.
- These instruments contain Nickel Titanium alloy. Do not use for individuals with known allergic sensitivity to these metals.
- For optimal cutting efficiently and safety do not use these instruments more than three times.

Precautions:

As with all products, you must exercise caution until you become proficient in its use. Length determination is imperative to ensure proper instrumentation using any rotary, reciprocating or hand instrument. The use of radiographs and/or apex locator are two acceptable methods of length determination. While we have implemented safeguards against possible misuse, there are several important points to remember:

- SafeSiders HF instruments are designed to dramatically increase cyclic fatigue resistance and improve torque strength. This design includes a helical flat which revolves in the opposite direction of the cutting flutes. This flat not only aids in the cutting process but assists in the removal of dentin coronally. This is not a manufacturing defect. Examine the instrument for any physical deformities before and during treatment. If any deformities are observed, instruments should be discarded.
- These instruments are to be used no more than three times. They can become less efficient after three uses causing undue stress to the instrument. This can lead to instrument separation. When navigating extremely curved canals, reduce instrument usage to one time.
- All instruments should be used with a gentle slow in-and-out pecking motion (2-3mm amplitude strokes). Work in the canal for the shortest possible time, and never use without irrigant in the canal.
- Always utilize minimal apical pressure. Don't force the instruments down the canal.
- These instruments are designed to easily go apically with minimal pressure.
- Operate the Endo-CSV™ handpiece at the 310 RPM (See Endo-CSV Instruction Manual).
- Operate the Endo-CSV™ handpiece at the #4 Torque setting (See Endo-CSV™ Instruction Manual). Straight-line access is a prerequisite for proper endodontic treatment. SafeSiders® HF instruments are no exception.
- Clean flutes after each 3 to 4 strokes during instrumentation.
- Frequently irrigate and lubricate the canal throughout the procedure.
- Exercise caution in the apical area and around significant curvatures.

Step-By-Step Instructions

For optimal results EDS recommends the use of SafeSiders® HF instruments in the Endo-CSV™ instrument in the "Auto-Reverse Forward" mode at 310 RPM set at 310 and the Torque settings set in the #4 position. See Endo-CSV instructions for complete details.

Caution, Consult Instructions for Use
NON-STERILE

Sterilization

- Instruments must be sterilized before use.
- Scrub the instruments with soap and warm water
- Rinse thoroughly with distilled or deionized water
- Allow to air dry
- Place the instruments in a sterilization pouch
- Use fresh distilled or deionized water
- Steam sterilize per gravity procedures below:

- Temperature: 270° F/132° C.
- Cycle Time: 15 minutes
- Minimum Dry Time: 30 minutes

- Inspect instruments prior to use for any abnormalities or deformities.
- Recommended Instruments Disposal: Place used instruments in Biohazard Sharps container.

In using the three instruments sequentially, negotiate to the apex starting with the Glidepath instrument (20/02), Shaper (30/02), Finisher (30/04). The progression of instrumentation and examination continues until you take the finishing instrument to working length.

Instrumentation Sequence

The basic SafeSiders HF System sequence will be the same, regardless of the canal size.

- Create straight-line access to canal orifice.
- Always irrigate. Establish working length by negotiating all root canals to their terminus with stainless steel hand files or reamers. EDS recommends the use of SafeSiders® hand instruments and an apex locator to achieve proper working length.
- Obtain the working length with a hand file preferably using an apex locator.
- Use the SafeSiders® HF Glide Path (Single Yellow Band) instrument in a gentle apical pecking motion, with short 2 – 3 mm amplitude strokes. Always utilize minimal apical pressure. Do not force the instruments down the canal. After 3 to 4 strokes, remove the instrument, inspect instruments for distortion, clean flutes and repeat irrigation until reaching the full working length.
- As entering the handpiece at 310 RPM (revolutions per minute).
- Finisher in the same manner as described above. It may be necessary to reconfirm patency and working length with a hand file (ISO 15), especially in more curved canals. After reaching the working length, each instrument can be used in a brushing motion away from external root concavities, to enlarge the root canal opening as well as buccally and lingually to ensure complete removal of pulp tissue in three dimensions.
- Finishing instruments should follow canal passively to working length, then be withdrawn.
- Confirm the finishing instrument reaches the apical measurement. The precision gutta percha point is now fit. (Note: Verify the gutta percha point fits to length).

Important Tips

- Always establish a reproducible glide path in each canal with a hand file or reamer (ISO 10 or higher) prior to rotary instrumentation.
- Irrigate, recapitulate and re-irrigate after using each instrument.
- Always confirm working length.
- Take instruments to resistance – never force instruments.
- Inspect instruments for distortion, wear, or bending.
- SafeSiders® HF instruments are to be used no more than three times.

Obturation

The SafeSiders® HF instrumentation system works with any obturation system (vertical, lateral condensation, master point, thermoplastic). EDS recommends the use of EZ-Fill® (EDS Cat. No. 1600-00) or EZ-Fill® Xpress (EDS Cat. No. 1620-00) root canal cement in conjunction with SafeSiders® HF precision gutta percha points. Pre-fit the SafeSiders HF 30/04 master gutta percha point (EDS Cat. No. 5104-G30) to ensure the point fits to length. In some cases you may need to select the SafeSiders HF 25/04 (EDS Cat. No. 5104-G25) point if working length is not reached with the 30/04.

NOTE: For more information, videos, and FAQ’s on this and all EDS products and techniques, please visit www.edsdental.com

For Use with Rotary Handpieces other than the Endo-CSV:

Indications for use: For the removal of dentin and shaping of the root canal.

Contraindications: None Known.

Warnings:

- Instruments are packaged non-sterile, and must be sterilized before use.
- These instruments contain Nickel Titanium alloy. Do not use for individuals with known allergic sensitivity to these metals.
- As endodontic rotary handpieces vary greatly, EDS recommends that the SafeSiders HF instruments are only used one time in endodontic handpieces other than the Endo-CSV™.

Precautions:

As with all products, you must exercise caution until you become proficient in its use. Length determination is imperative to ensure proper instrumentation using any rotary, reciprocating or hand instrument. The use of radiographs and/or apex locator are two acceptable methods of length determination. While we have implemented safeguards against possible misuse, there are several important points to remember:

- SafeSiders HF instruments are designed to dramatically increase cyclic fatigue resistance and improve torque strength. This design includes a helical flat which revolves in the opposite direction of the cutting flutes. This flat not only aids in the cutting process but assists in the removal of dentin coronally. This is not a manufacturing defect. Examine the instrument for any physical deformities before and during treatment. If any deformities are observed, instruments should be discarded.
- As entering the rotary handpiece vary greatly, EDS recommends that the SafeSiders HF instruments are only used one time in endodontic handpieces other than the Endo-CSV™.
- All instruments should be used with a gentle slow in-and-out pecking motion (2-3mm amplitude strokes). Work in the canal for the shortest possible time, and never use without irrigant in the canal.
- These instruments are designed to easily go apically with minimal pressure.
- Always utilize minimal apical pressure. Don't force the instruments down the canal.
- A slow-speed handpiece is required for rotary file use.
- Operate the handpiece at 310 RPM (revolutions per minute).
- Straight-line access is a prerequisite for proper endodontic treatment. SafeSiders HF instruments are no exception.
- Clean flutes after each 3 to 4 strokes during instrumentation.
- Frequently irrigate and lubricate the canal throughout the procedure.
- Exercise caution in the apical area and around significant curvatures.

Step-By-Step Instructions

SafeSiders HF instruments are used in a controlled, slow-speed rotary handpiece (set at 310 RPM), thus providing a continuous tapered preparation from orifice to apex but mirrors the original canal anatomy. A torque control motor is recommended.

Sterilization

Instruments must be sterilized before use.

- Scrub the instruments with soap and warm water
- Rinse thoroughly with distilled or deionized water
- Allow to air dry
- Place the instruments in a sterilization pouch
- Use fresh distilled or deionized water
- Steam sterilize per gravity procedures below:

- Temperature: 270° F/132° C.
- Cycle Time: 15 minutes
- Minimum Dry Time: 30 minutes

- Inspect instruments prior to use for any abnormalities or deformities.
- Recommended Instruments Disposal: Place used instruments in Biohazard Sharps container.

In using the three instruments sequentially, negotiate to the apex starting with the Glidepath instrument (20/02), Shaper (30/02), Finisher (30/04). The progression of instrumentation and examination continues until you take the finishing instrument to working length.

Instrumentation Sequence

The basic SafeSiders HF System sequence will be the same, regardless of the canal size.

- Create straight-line access to canal orifice.
- Always irrigate. Establish working length by negotiating all root canals to their terminus with stainless steel hand files or reamers. EDS recommends the use of SafeSiders hand instruments and an apex locator to achieve proper working length.
- Obtain the working length with a hand file (ISO 10 or higher) preferably using an apex locator.
- Use the SafeSiders HF Glide Path (Single Yellow Band) instrument in a gentle apical pecking motion, with short 2 – 3 mm amplitude strokes. Always utilize minimal apical pressure. Do not force the instruments down the canal. After 3 to 4 strokes, remove the instrument, inspect instruments for distortion, clean flutes and repeat irrigation until reaching the full working length.
- Use the SafeSiders HF Shaper (Single Blue Band) and SafeSiders HF (Double Blue Band) Finisher in the same manner as described above. It may be necessary to reconfirm patency and working length with a hand file (ISO 15), especially in more curved canals. After reaching the working length, each instrument can be used in a brushing motion away from external root concavities, to enlarge the root canal opening as well as buccally and lingually to ensure complete removal of pulp tissue in three dimensions.
- Finishing instruments should follow canal passively to working length, then be withdrawn.
- Confirm the finishing instrument reaches the apical measurement. The precision gutta percha point is now fit. (Note: Verify the gutta percha point fits to length).

Important Tips

- Always establish a reproducible glide path in each canal with a hand file or reamer (ISO 10 or higher) prior to rotary instrumentation.
- Irrigate, recapitulate and re-irrigate after using each instrument.
- Always confirm working length.
- Take instruments to resistance – never force instruments.
- Inspect instruments for distortion, wear, or bending.
- SafeSiders HF instruments are only to be used once in rotation.

Obturation

The SafeSiders® HF instrumentation system works with any obturation system (vertical, lateral condensation, master point, thermoplastic). EDS recommends the use of EZ-Fill® (EDS Cat. No. 1600-00) or EZ-Fill® Xpress (EDS Cat. No. 1620-00) root canal cement in conjunction with SafeSiders® HF precision gutta percha points. Pre-fit the SafeSiders HF 30/04 master gutta percha point (EDS Cat. No. 5104-G30) to ensure the point fits to length. In some cases you may need to select the SafeSiders HF 25/04 (EDS Cat. No. 5104-G25) point if working length is not reached with the 30/04.

NOTE: For more information, videos, and FAQ’s on this and all EDS products and techniques, please visit www.edsdental.com

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89 Leuning Street, S. Hackensack, NJ 07606
www.edsdental.com

U.S. and other foreign patents pending.

  European Healthcare & Devices Ltd., Stratton House, Bishopstown Rd, Cork, T12 Y9TC, Ireland. Contact: TQMUK@aol.com

Caution: Federal law restricts this device to sale by or on the order of a licensed dentist.

Important Sterilization Procedures: **Endodontic Instruments, Posts, Drills and** **Wrenches are non-sterile.**

Prior to use, bag materials in an autoclave safe pouch and sterilize with a gravity displacement autoclave for 15 minutes at 132°C with a minimum drying time of 30 minutes. When reprocessing instruments, drills, wrenches, or taps, remove debris and dry before sterilization. Remove debris/soil using the EDS recommended cleaning protocol. Disinfect using a thermal disinfectant unit (washer-disinfectant) at 90°C for 1 minute.

For complete instructions see website.

89 Leuning Street, S. Hackensack, NJ 07606
Tel: 1-800-22-FLEXI, 201-487-9090 • Fax: 201-487-5120
Web: www.edsdental.com E-Mail: info@edsdental.com

Non-Sterile

(OVER)