

CEMPOWER™

Dual-Cured, Self-Etch, Self-Adhesive Resin Cement

CEMPOWER is a self-etch, self-adhesive, dual cured resin cement. **This user-friendly cement is ideal for restorations made of zirconium oxide, lithium disilicate, other metal-ceramics, and for the cementation of high-strength restoratives on implant abutments.**



Includes: 8 gm Auto Mix Syringe,
5 Straight Mixing Tips, 5 Angled Intra canal tips.
See back for details.

Available in A2 and Translucent Shades

Feature


- New adhesion promoter technology ensures high retention on multiple surfaces
- Superior Hydrolytic Stability
- Numerous Restorative uses
- Color Fast Proprietary Initiator System
- Biocompatible, Bis-GMA Free, Patient Friendly Self Adhesive Chemistry
- “Tack and Wave” Chemistry
- Minimal Oxygen Inhibition
- Time-saving Technique
- Stable Material Properties
Mixes quickly, Dual cured

Benefits

- Flexible use, high retention of zirconia, lithium disilicate, enamel, dentin, porcelain fused to metal (PFM), gold, etc.
- Ensures long-lasting highly retentive restoration
- Use for inlays, onlays, crowns, bridges, posts and veneer cementation
- Color stability to extend the life of your aesthetic restoration
- Foremost patient wellbeing with virtually no post-op sensitivity
- Enables restorations to be initially light-cured into place to establish a “semi-gel” state for spot-tacking into place and easy cleaning of the margins
- Ensures curing stability of material
- No etching, priming or bonding - Saves you time and money
- No refrigeration needed, long shelf life
- Simple, single step cementation


THE POWER OF CEMPOWER

CEMPOWER Cement




Bond Strength To Dentin

Bond Strength To Enamel



Bond Strength To Zirconium



ORDERING INFORMATION

CEMPOWER Cement:

- 930-00: A2 Shade: 8gm Auto Mix Syringe, 5 Straight Mixing Tips, 5 Angled Intra canal tips, 1 mixing pad.
- 940-00: Translucent Shade: 8gm Auto Mix Syringe, 5 Straight Mixing Tips, 5 Angled Intra canal tips, 1 mixing pad.

Mixing Tip Refills

- 870-20: 20 - Angled Intra canal tips.
- 1628-20: 20 - Straight Mixing Tips

For our complete line of patented dental products please visit our website at: www.edsdental.com *Data On File